

B'SYS GmbH

CHO $\alpha_1\beta_2\gamma_2$ GABA_A Cell Line

Specification Sheet

© B'SYS GmbH

TABLE OF CONTENTS

1	BACKGROUND	3
1.1	THE PHARMACOLOGICAL DISTINCTION OF GABA _A RECEPTOR SUBTYPES	3
1.2	B'SYS' CHO $\alpha_1\beta_2\gamma_2$ CELLS	3
2	PRODUCT SHIPMENT	3
2.1	PRODUCT FORMAT	3
2.2	MYCOPLASMA CERTIFICATE	3
3	VALIDATION OF CHO $\alpha_1\beta_2\gamma_2$ GABA_A CELLS	4
3.1	ELECTROPHYSIOLOGY	4
3.2	GABA AS AGONIST	4
3.3	POSITIVE ALLOSTERIC AGONISTS (PAM)	4
3.4	ANTAGONIST	5
3.5	PATCH-CLAMP SUCCESS RATES	6
4	CELL CULTURE CONDITIONS	6
4.1	GENERAL	6
4.2	RECOMMENDED COMPLETE MEDIUM	6
4.3	ANTIBIOTICS	6
4.4	THAWING CELLS	6
4.5	SPLITTING CELLS	7
4.6	FREEZING MEDIUM	7
4.7	FREEZING CELLS	7
5	SEQUENCE	8
5.1	GABA α_1 , NP_000797.2	8
5.2	GABA β_2 , NP_000804.1	8
5.3	GABA γ_2 , NP_000807.2	8
6	CONTACT INFORMATION	8
6.1	CONTACT ADDRESS FOR TECHNICAL SUPPORT & ORDERING INFORMATION	8

1 BACKGROUND

1.1 The Pharmacological Distinction of GABA_A Receptor Subtypes

Using genetically modified (knock-in) mice it has been demonstrated that GABA_A receptors containing an α_1 subunit mediate the sedative/muscle relaxant effects of benzodiazepines, whereas β_2 and/or β_3 subunit containing receptors mediate the anxiolytic and anticonvulsant effects. GABA_A α_5 receptors have a relatively restricted distribution being primarily expressed in the hippocampus, a region of the brain associated with learning and memory, and although α_5 receptors account for less than 5% of the total GABA_A receptor population in the brain, in the hippocampus they represent 20% of all GABA_A receptors, thereby implicating this GABA_A receptor subtype in learning and memory processes. Thus, the pharmacological distinction of GABA_A receptor isoforms serves now as a promising basis for the development of new agents effective at only restricted brain regions and thus exhibiting unique and specific physiological effects.

1.2 B'SYS' CHO $\alpha_1\beta_2\gamma_2$ Cells

B'SYS has designed a CHO $\alpha_1\beta_2\gamma_2$ GABA_A receptor cell line with constitutive expression of human α_1 subunit together with the β_2 and γ_2 subunits. The human GABA_A receptor cDNAs were cloned and transfected into CHO cells and then the functional properties of the GABA_A receptors validated by means of manual and automated (Q-Patch) patch-clamp technique. Results are outlined in section 3.

2 PRODUCT SHIPMENT

2.1 Product Format

CHO cells stably transfected with recombinant $\alpha_1\beta_2\gamma_2$ human GABA_A receptors:

- 0.5 mL aliquots of frozen cells at approximately 3.0 E+06 cells/mL
- Cells are frozen in complete medium with 10% DMSO

2.2 Mycoplasma Certificate

B'SYS periodically tests cells for presence of mycoplasma by means of highly sensitive PCR based assays. All delivered cells are free of mycoplasma.

3 VALIDATION OF CHO $\alpha_1\beta_2\gamma_2$ GABA_A CELLS

3.1 Electrophysiology

GABA_A currents were measured by means of the patch-clamp technique in the whole-cell configuration. The bath solution contained (in mM) NaCl 137, KCl 4, CaCl₂ 1.8, MgCl₂ 1, D-glucose 10, HEPES 10, pH (NaOH) 7.40. The pipette solution consisted of (in mM) KCl 130, MgCl₂ 1, MgATP 5, HEPES 10, EGTA 5, pH (KOH) 7.20. After formation of a GΩ seal between the patch electrodes and individual $\alpha_1\beta_2\gamma_2$ GABA_A stably transfected CHO cells, the cell membrane across the pipette tip was ruptured to assure electrical access to the cell interior. All solutions applied to cells were maintained at room temperature. The holding potential was -80 mV. As soon as a stable seal could be established inward chloride currents were measured upon application of agonist concentrations (GABA) to patch-clamped cell.

3.2 GABA as Agonist

The concentration dependence for GABA was tested using automated patch-clamping (Q-Patch). Concentrations between 0.01 μM and 100 μM were applied. Increasing concentrations were tested every 30 s. GABA was washed off after each GABA application. The EC₅₀ was determined to be 5.92 μM (Hill coefficient 2.4).




Fig.1: Activation of GABA_A receptor currents upon brief applications of increasing concentrations of GABA to patch-clamped cell. The cell membrane was held at -80 mV. No inward chloride currents were recorded in untransfected cells (data not shown).

3.3 Positive allosteric agonists (PAM)

As positive allosteric agonists Diazepam and Propofol were tested in concentrations between 0.01 μM and 30 μM in the presence of 2.0 μM GABA:

Diazepam: EC₅₀: 0.11 μM, Hill: 1.29, a_{max}: 152%

Propofol: EC₅₀: 3.59 μM, Hill: 1.31, a_{max}: 763%




Fig.2: Effects of positive allosteric modulators (in the presence of 2.0 μM GABA): Top Diazepam, Bottom Propofol

3.4 Antagonist

As antagonist Bicuculline was tested at concentration between 0.01 μM and 10 μM in the presence of 5.0 μM GABA. The IC_{50} was determined to be: 0.48 μM (Hill coefficient: 0.92).




Fig.3: Effects of the antagonist Bicuculline (in the presence of 5.0 μM GABA)

3.5 Patch-clamp Success Rates

The patch-clamp properties of the CHO $\alpha_1\beta_2\gamma_2$ GABA_A receptor cell line were elucidated at typical working passage numbers (passage 2-18). A total of 63 cells were analyzed. Success for establishment of on-cell configuration was defined as follows: > 1 G Ω . The whole-cell configuration was not accepted if the membrane resistance was below 500 M Ω . A successful recording had to be free of rundown effects and variations in series resistance.

- On-cell successful: **95%** (n=80)
- Whole-cell successful: **91%** (n=80)
- Recording (15 min) successful: **78%** (n=80)

4 CELL CULTURE CONDITIONS

4.1 General

CHO $\alpha_1\beta_2\gamma_2$ GABA_A receptor cells are incubated at 37°C in a humidified atmosphere with 5% CO₂ (rel. humidity > 95%). The cells are continuously maintained and passaged in sterile culture flasks containing F12 (HAM) medium supplemented with 10% fetal bovine serum, 1.0% Penicillin/Streptomycin solution and Hygromycin 200 μ g/mL, Puromycin 5 μ g/mL, Zeocin 100 μ g/mL. The CHO $\alpha_1\beta_2\gamma_2$ GABA_A receptor cells are passaged at a confluence of about 50-80%. For electrophysiological measurements the cells are seeded onto e.g. 35 mm sterile culture dishes containing complete medium.

- All solutions and equipment coming in contact with the cells must be sterile.
- Use proper sterile technique and work in a laminar flow hood.
- Be sure to have frozen cell stocks at hand before starting experiments.
- Cells should be split every 4-5 days at 50% - 80% confluency at 1:3 to 1:5 ratio.

4.2 Recommended Complete Medium

- DMEM/F12 with L-Glutamine or GlutaMAX I
- 10% FBS
- 1.0% Penicillin/Streptomycin

4.3 Antibiotics

- CHO $\alpha_1\beta_2\gamma_2$ •GABA_A receptor clones were selected under Hygromycin 500 μ g/mL, Puromycin 5 μ g/mL, Zeocin 100 μ g/mL antibiotic pressure.
- To cultivate $\alpha_1\beta_2\gamma_2$ GABA_A receptor cells, a reduced antibiotic pressure (Hygromycin 200 μ g/mL, Puromycin 5 μ g/mL, Zeocin 100 μ g/mL) should be used.

Remark: The permanent application of antibiotic pressure has no effect on current density.

4.4 Thawing Cells

- Remove vial of cells from liquid nitrogen and thaw quickly at 37°C.
- Decontaminate outside of vial with 70% ethanol.
- Transfer cells to a T-25 culture flask containing 5 mL complete medium with Hygromycin 200 μ g/mL, Puromycin 5 μ g/mL, Zeocin 100 μ g/mL.
- Incubate cells and check them daily until 50% - 80% confluency is reached.

4.5 Splitting Cells

- When cells are 50% - 80% confluent remove complete medium.
- Wash cells once with 1x Trypsin/EDTA.
- Remove Trypsin/EDTA quickly and incubate cells for 3-5 min at room temperature.
- Detach cells, add complete medium and pipet up and down to break clumps of cells.
- Passage cells into new flask with complete medium and antibiotics at 1:3 to 1:5 ratio.
- Use remaining suspension for counting the cells.

4.6 Freezing Medium

- Mix 0.9 mL fresh complete medium and 0.1 mL DMSO for every 1 mL freezing medium.
- Sterilize freezing medium by means of appropriate micro filter (0.1 μm – 0.2 μm).

4.7 Freezing Cells

- Prepare fresh freezing medium and keep it on ice.
- Cells should have 80% - 90% confluency prior to freezing.
- Remove the complete medium.
- Wash cells once with 1x Trypsin/EDTA.
- Remove Trypsin/EDTA quickly and incubate cells for 3-5 min at room temperature.
- Detach cells, add complete medium and pipet up and down to break clumps of cells.
- Pellet cells with centrifuge and carefully aspirate off medium.
- Resuspend cells at a density of approximately 3.0×10^6 cells per mL with fresh freezing medium.
- Aliquot 0.5 mL of cell suspension into each cryovial.
- Overnight incubate cells in a styropor box at -80°C .
- The next morning transfer cryovial in liquid nitrogen tank for long-term storage.

5 SEQUENCE

Some subunits were codon optimized.

5.1 GABA α_1 , NP_000797.2

MRKSPGLSDCLWAWI LLLSTLTGRSYGQPSLQDELKDNTTTFTRILDRLLDGYDNRLRPLGLGERVTEVKTDIFVTSFGPV
SDHDMEYTI DVFFRQSWKDERLKFKGPM TVLRLNNLMASKI WTPDTFFHNGKKSVAHNMTMPNKLRLRITEDGTLLYTMRL
TVRAECPMHLEDFPMDAHACPLKFGSYAYTRAEVVYEW TREPARSVVVAEDGSRLNQYDLLGQTVDSGIVQSSTGEYVVM
TTHFHLKRKIGYFVIQTYLPCIMTVILSQVSFWLNRESVPARTVFGVTTVLTMTTLSISARNSLPKVAYATAMDWFI AVC
YAFVFSALIEFATVNYFTKRGYAWDGKSVVPEKPKKVKDPLIKKNNTYAPTATSYPNLRAGDPGLATIAKSATIEPKEV
KPETKPEPKKTFNSVSKIDRLSR IAFPLLFGIFNLVYWATYLNREPQLKAPTPHQ*

5.2 GABA β_2 , NP_000804.1

MWRVRKRGYFGIWSFPLIIAAVCAQSVNDPSNMSLVKETVDRLLKGYDIRLRPDPFGGPPVAVGMNIDIASIDMVSEVNMD
YTLTMYFQQAWRDKRLSYNVIPLNLTLDNRVADQLWVPDQYFLNDKKS FVHGVTVKNRMIRLHPDGTVLYGLRITTTAAC
MMDLRRYPLDEQNCTLEIESYGYTTDDIEFYWRGDDNAVTVGVTKIELPQFSIVDYKLITKKVVFSTGSYPRLSLSFKLKR
NIGYFILQTYMPSILITILSWVSFWINYDASAARVALGITTVLTMTTINTHLRETLPKIPYVKAIDMYLMGCFVVFVFMAL
LEYALVNYIFFGRGPQRQKKA AEKAASANNEKMRLDV NKMDPHENILLSTLEIKNEMATSEAVMGLGDPRSTMLAYDASS
IQYRKAGLPRHSFGRNALERHVAQKKSRLRRRASQLKITIPDLTDVNAIDRWSRIFFPVVFSFFNIVYWLYYVN*

5.3 GABA γ_2 , NP_000807.2

MSSPNIWSTGSSVYSTPVFSQKMTVWILLLLSLYPGFSTQKSDDDYEDYASNKTWVLT PKVPEGDVTVILNNLLEGYDNK
LRPDIGVKPTLIHTDMYVNSIGPVNAINMEYTI DIFFAQTWYDRRLKFNSTIKVLRLNSNMVGIWIPDTFFRNSKKADA
HWITTPNRMLRIWNDGRVLYTLRLTIDAECQLQLHNFPMDEHSCPLEFSSYGYPREEIVYQWKRSSVEVGDTRSWRLYQF
SFVGLRNTTEVVKTTSGDYVMSVYFDLSRRMGYFTIQTYIPCTLIVVLSWVSFWINKDAVPARTSLGITTVLTMTTLLST
IARKSLPKVSYVTAMDLFVSVCFIFVFSALVEYGT LHYFVSNRKP SKDKDKKKKNPAPTIDIRPRSATI QMNNATHLQER
DEEYGYECLDGKDCASFFCCFEDCRTGAWRHGRIHIRIAKMDSYARIFFPTAFCLFNLVYVWSYLYL*

6 CONTACT INFORMATION

6.1 Contact Address for Technical Support & Ordering Information

- B'SYS GmbH
Technology Center Witterswil
Benkenstrasse 254
4108 Witterswil
Switzerland

Tel: +41 61 721 77 44

Fax: +41 61 721 77 41

Email: info@bsys.ch

Web: www.bsys.ch